

A novel and noninvasive FIT, DNA, mRNA and machine learning / AI-based early Colorectal Cancer and Advanced Adenoma detection test

First interim data review of the international COLOFUTURE study

L. Krammes*, A. Jüngling*, H. Mahmood*, F. Frondorf*, J. Luther*, V. Königs*, C. Schölz*, P. Becker*, S. Maharjan*, A. Sever*, S. Garapati*, A. Balasubramaniam*, M. Knabe*, R. Lovlie*, **M. Eidens***, M. Dollinger**

* Mainz Biomed Germany GmbH, Mainz, Germany
** MD, Professor of Gastroenterology, Klinikum Landshut Anstalt des öffentlichen Rechts der Stadt Landshut, Germany



INTRODUCTION

Colorectal cancer (CRC) ranks as the **second leading cause** of global cancer mortality, with incidence particularly rising among younger demographics. Early detection of polyps and lesions that can progress to CRC at the **advanced adenoma (AA) stage** or earlier is crucial as it prevents them from progressing further and leads to better treatment outcomes. As colonoscopy participation rates remain poor, reliable noninvasive screening methods like Fecal Immunochemical Test (FIT) or novel alternatives are needed. Here, we describe a **new approach** applying state-of-the-art PCR based **DNA and mRNA-based screening** in combination with FIT and a ML/AI driven algorithm, leading to a highly sensitive and specific clinical test for CRC and AA.


OBJECTIVES

- Interim analysis of the clinical performance of a novel noninvasive FIT, DNA and mRNA-based CRC screening test using a proprietary reagent stabilization formulation, small stool sample volumes, state-of-the-art machine learning (ML) and an artificial intelligence (AI) driven algorithm.
- Demonstrate high-confidence clinical classification based on minimal sample input.
- Discrimination of CRC and its precursor lesions (including AAs) from normal control samples with high sensitivity and specificity.

METHODS

COLOFUTURE Trial Overview

- Prospective, multi-center study across Germany and Norway.
- Study participants were 40-85 years old, referred for colonoscopy (screening or diagnostic), or had a diagnosis of colorectal adenocarcinoma with no previous treatment.



Laboratory Process Overview




Figure 1: COLOFUTURE laboratory process overview.

- Occult blood was determined by hemoglobin quantification.
- The proprietary DNA/mRNA stool collection tube and sample underwent an automated silica-bead based extraction method followed by (RT)-qPCR for novel mRNA and DNA biomarkers.
- The mRNA biomarkers were analyzed by absolute and relative quantification approaches. For relative quantification, normalization by a housekeeping gene (GAPDH) was used.
- The artificial intelligence/ machine learning algorithm was developed in partnership with Liquid Biosciences of Aliso Viejo, California.

COLOFUTURE Patient Age and Diagnosis


| | Mean age | Control | Non-AA | AA | CRC | All groups |
|-------------|----------|---------|--------|----|-----|------------|
| Male | 62.2 | 56 | 6 | 37 | 23 | 122 |
| Female | 62.4 | 61 | 4 | 18 | 15 | 98 |
| All genders | 62.3 | 117 | 10 | 55 | 38 | 220 |

RESULTS

Genetic Biomarkers Analyzed

| Biomarker | Name | Function | Target |
|---------------|---|--|--------|
| BRAF, somatic | Serine/ threonine kinase | MAP kinase/ ERK signaling pathway | DNA |
| KRAS, somatic | Kirsten ras oncogene homolog | GTPase, regulation of cell proliferation | DNA |
| B2M | Beta-2-microglobulin | Antigen presentation in immune system | mRNA |
| CEACAM5 | Carcinoembryonic antigen-related cell adhesion molecule 5 | Cell adhesion, Intracellular signaling | mRNA |
| ITGA6 | Integrin alpha-6 heavy chain | Cell surface adhesion and signaling | mRNA |
| MACC1 | Metastasis-associated in colon cancer protein 1 | Transcription activator | mRNA |
| PTGS2 | Prostaglandin G/H synthase 2 | Production of inflammatory prostaglandins | mRNA |
| S100A4 | S100 calcium binding protein A4 | Regulation of cellular processes such as cell cycle progression and differentiation | mRNA |
| GAPDH | glyceraldehyde-3-phosphate dehydrogenase | Moonlighting protein based on its ability to perform mechanically distinct functions | mRNA |

Distribution of Pathological Results for CRC and AA



Clinical Performance

Performance of advanced nucleic acid screening for AA and CRC combined with classical FIT. Performance is displayed as sensitivity and specificity for CRC, AA and a combined group of diseased subjects with two-sided 95% Clopper-Pearson confidence intervals (CI).

| Category | DNA+mRNA+FIT | | FIT only* | |
|--|----------------------------|----------------------------|----------------------------|----------------------------|
| | Sensitivity (%) | Specificity (%) | Sensitivity (%) | Specificity (%) |
| CRC vs. Normal + non-AA (95% confidence interval) | 94.4 (81.3-99.3) | 97.5 (92.9-99.5) | 73.7 (56.9-86.6) | 93.7 (88.0-97.2) |
| AA vs. Normal + non-AA (95% confidence interval) | 80.0 (66.3-90.0) | 95.2 (89.8-98.2) | 29.1 (17.6-42.9) | 93.7 (88.0-97.2) |
| CRC + AA vs. Normal + non-AA (95% confidence interval) | 87.1 (78.5-93.2) | 92.9 (87.0-96.7) | 47.3 (36.9-57.9) | 93.7 (88.0-97.2) |

*calculated on all available samples

CRC Stage Sensitivities

| Stage | Pathology staging* | Detected | Sensitivity (%) |
|---------|--------------------|----------|-----------------|
| Stage 0 | 1/34 | 0 | 0 |
| Stage 1 | 4/34 | 4 | 100 |
| Stage 2 | 16/34 | 15 | 94 |
| Stage 3 | 11/34 | 11 | 100 |
| Stage 4 | 2/34 | 2 | 100 |

* 4 CRC samples were not staged at the time of the interim analysis

Sensitivity for high grade dysplasia: 75% (9/12)

CONCLUSION

Applying this novel approach, **sensitivity** for detection of **CRC** was found to be **94.4%** while **precursor lesions including AAs** were detected with **80% sensitivity**. Specificity (negative colonoscopies plus non-AAs) was calculated to be 97.5% and 95.2% for CRC and AA, respectively. This innovative setup represents the inaugural instance of a multimodal analysis involving DNA, mRNA expression, FIT results and an AI/ML developed algorithm. Results indicate a substantial and meaningful enhancement in the effectiveness of non-invasive CRC screening, particularly for the detection of AA, where increased sensitivity is urgently needed to significantly reduce CRC incidence and mortality.

References:

- Herring E, Kanaoka S, Tremblay E, Beaulieu JF. A Stool Multitarget mRNA Assay for the Detection of Colorectal Neoplasms. *Methods Mol Biol*. 2018;1765:217-227. doi: 10.1007/978-1-4939-7765-9_14. PMID: 29589311.
- Herring E, Tremblay E, McFadden N, Kanaoka S, Beaulieu JF. Multitarget Stool mRNA Test for Detecting Colorectal Cancer Lesions Including Advanced Adenomas. *Cancers (Basel)*. 2021 Mar 11;13(6):1228. doi: 10.3390/cancers13061228. PMID: 33799738; PMCID: PMC7998137.
- Gies A, Cuk K, Schrotz-King P, Brenner H. Direct comparison of ten quantitative fecal immunochemical tests for hemoglobin stability in colorectal cancer screening. *Clin Transl Gastroenterol*. 2018 Jul 6;9(7):168. doi: 10.1038/s41424-018-0035-2. PMID: 29976921; PMCID: PMC6033915.
- Gies A, Niedermaier T, Gruner LF, Heisser T, Schrotz-King P, Brenner H. Fecal Immunochemical Tests Detect Screening Participants with Multiple Advanced Adenomas Better than T1 Colorectal Cancers. *Cancers (Basel)*. 2021 Feb 5;13(4):644. doi: 10.3390/cancers13040644. PMID: 33562775; PMCID: PMC7914536.
- Dekker E, Tanis PJ, Vieugels JLA, Kasi PM, Wallace MB. Colorectal cancer. *Lancet*. 2019 Oct 19;394(10207):1467-1480. doi: 10.1016/S0140-6736(19)32319-0. PMID: 31631858.
- Ladabaum U, Church TR, Feng Z, Ransohoff DF, Schoen RE. Counting Advanced Precancerous Lesions as True Positives When Determining Colorectal Cancer Screening Test Specificity. *J Natl Cancer Inst*. 2022 Jul 11;114(7):1040-1043. doi: 10.1093/jnci/djac027. PMID: 35134969; PMCID: PMC9275773.

Contact details:

Dr. Moritz Eidens
Mainz Biomed N.V.
Sirius Gutenberg Park
Robert-Koch-Straße 50
55129 Mainz, Germany



Email Address: Moritz.eidens@mainzbiomed.com